


# BACTERIOLOGICAL INCUBATOR



## FEATURES

### Optional / Additional with Instruments (if Required)

- Trays shift as per requirement.
- Microprocessor based PID temperature controller/
- Microprocessor Based PID Controller alongwith PC communication thru built-in RS-485 interface and printing option through software/
- Touch Screen Display (3.5") system & Pen Drive System alongwith data logging & for saving of data and communication to PC.

## SPEC. / MODEL

MODEL NO.	SIZE INSIDE CHAMBER ( L X W X H )	LTRS. APPROX.
	DIA X DEPTH MM	
BioBAC-028	12" X 12" X 12" (300 X 300 X 300 MM)	28 LTRS.
BioBAC-045	14" X 14" X 14" (350 X 350 X 350 MM)	45 LTRS.
BioBAC-095	18" X 18" X 18" (450 X 450 X 450 MM)	95 LTRS.
BioBAC-0128	18" X 18" X 24" (450 X 450 X 600 MM)	128 LTRS.
BioBAC-0228	24" X 24" X 24" (600 X 600 X 600 MM)	228 LTRS.
BioBAC-0256	24" X 18" X 36" (600 X 450 X 900 MM)	256 LTRS.
BioBAC-0342	24" X 24" X 36" (600 X 600 X 900 MM)	342 LTRS.

Instrument Manufactured and supplied as per required dimension.

**Delhi Office:**  
S-10, Pankaj Plaza,  
Plot No.5, Pocket-5,  
Sector 12, Dwarka,  
New Delhi-110078  
(INDIA)

**Contacts:**  
sales@biogenixsystems.com  
biogenixsystems@gmail.com  
+91-9911717057, 9716064799  
9990251199, 011-43738899  
www.biogenixsystems.com

○ **MARKETED BY**

○ **BIOGENIX SYSTEMS**

○ (AN ISO 9001:2008 CERTIFIED CO.)

- Temp. Controlled by Digital Temp. controller.
- Inner Chamber double walled Stainless Steel Outer Body M.S. Duly epoxy powder coated/Stainless Steel.
- Stainless steel with perforated adjustable shelves.
- Specially designed to maintain 37°C
- Heating elements are placed in ribs at the bottom sides.
- Temp. Range +5°C to 70°C
- Temperature accuracy : ± 0.3 C at 37° C
- Temperature Increment ±0.1° C
- CFC Free Compressor of "Kirloskar",
- Full view inner glass door. Lock & Key arrangement,
- Caster wheel mounted cabinet
- The inner chamber is properly insulated with special PUF to prevent thermal losses.
- Two Nos of Air Circulating Fan for uniform temperature, alongwith LED light illumination system inside the chamber.
- Trays are arranged as per chamber dimension.

**WARRANTY: ONE YEAR ADDITIONAL YEARS AS REQUIRED**