

LABORATORY VERTICAL AUTOCLAVE

BRAND - BIOGENIX®

Product Description and Specification

For catering the variegated needs of our clients, we are deeply engrossed in offering a qualitative array of **Laboratory Vertical Autoclave**.

KEY FEATURES

- Stainless steel structure
- Hand wheel type / Screw type
- Door safety lock system
- Digital Display
- Automatic exhaust system for cold air and pressure Over-temperature and over-pressure auto-protection
- Alarm function for completion of process and low water level Self-inflating type seal
- Optional: Drying system



SPECIAL FEATURES

- Computer controlled auto recycle sterilization.
- LCD screen indicating working state,
- hand wheel type for a quick door opening.
- pressure safety lock.
- Temp. scope: 50-126 Deg. C
- Time scope: 0-99 hours.
- Mini printer process recorder can be attached - Optional.
- Power: 2500 Watt.
- Includes baskets.

SPECIFICATION

Model	BS-WAV-101	BSWAV-102	BSWAV-103	BSWAV-104
Capacity	35L	50 L	75 L	100 L
Temperature Range	105~134°C			
Sterilization Temperature	134°C			
Heat Average	≤±1°C			
Working Pressure	0.22 MPa			
Maximum Working Pressure	0.23 Mpa			
Timer Range	0~60 min			
Overall Dimension (LxWxH)	480x460x850 mm	520x520x980 mm	560x560x980 mm	590x590x1080 mm
Packing Size	570x550x970 mm	590x590x1110 mm	650x630x1150 mm	680x650x1220 mm
Weight	42 / 56 kg	50 / 68 kg	70 / 90 kg	85 / 105 kg
Power	2500 W	3000 W	4500 W	
Power Supply	AC220V 50Hz			

Size, Capacity & ordering Information

Model	Size / Capacity (Litres) Including Lid	
BSWVA-101	570x550x970 mm	35L
BSWVA-102	590x590x1110 mm	50 L
BSWVA-103	650x630x1150 mm	75 L
BSWVA-104	680x650x1220 mm	100 L



BIOGENIX SYSTEMS PRIVATE LIMITED

S-10, Pankaj Central Plaza, Pocket-5, Plot no-5,
Sector-12, Dwarka, New Delhi-110075, INDIA

Email: sales@biogenixsystems.com |

Phone: +91-11-43738899 | 9990251155-99 | 9911717057

Website: www.biogenixsystems.com

Dealer Details