

ORBITAL SHAKER

Brand- **BIOGENIX®**



Product Description

Our product is an accurately designed, microprocessor controlled with double Decker platform to save your valuable lab space. Three eccentric shaft balancing drive ensures shaking with uniform speed. Speed adjustment settings permit both gentle and vigorous shaking. Varieties of platforms are available for different glassware and vessels.



BSOT-601

600 SERIES ORBITAL SHAKER

Microprocessor based PID controller with timing function

Brushless DC motor ensures maintenance-free, reliable and quiet operation

Easy to clean anti-eroding stainless steel or alnico plate Swift start and stop system prevents liquid spillage

Audio alarm signals program completion

SPECIFICATIONS:

Model	BSOT-601	BSOT-602	BSOT-603	BSOT-604
Speed Range	50-250 rpm		40-300 rpm	
Temperature	No heating function			
Shaking Orbit	10 mm	20 mm		
Timing Range	1~9999 min			
Motion	Orbital			
Operation	Continuous/Timed			
Platform Dimension	280x270 mm	350x260 mm	400x340 mm	
Display	Analog			
Power	100 W	110 W	120 W	
Power Supply	220V, 50Hz			

Model	BSOT-605	BSOT-606	BSOT-607	BSOT-608
Speed Range	40-300 rpm	20-250 rpm		40-300 rpm
Temperature	No heating function	RT 5-100°C		RT+5-65°C
Shaking Orbit	20 mm			
Timing Range	1~9999 min			
Motion	Orbital			
Operation	Continuous/ Timed			
Platform Dimension	450x450 mm	350x260 mm	380x430 mm	400x340 mm
Display	Analog	LC D		
Power	130 W	300 W	400 W	450 W
Power Supply	220V, 50Hz			



BSOT-602



BSOT-603

Dealer Details



BIOGENIX SYSTEMS PRIVATE LIMITED

S-10, Pankaj Central Plaza, Pocket-5, Plot no-5,
Sector-12, Dwarka, New Delhi-110075, INDIA

Email: sales@biogenixsystems.com |

Phone: +91-11-43738899 | 9911717057

Website: www.biogenixsystems.com

info@biogenixsystems.com
www.biogenixsystems.com