

# Rotating Mixture

## Brand - BioGenix®

### PRODUCT DESCRIPTION

**BIOGENIX®** Rotating Mixture is a General purpose instrument which is designed for the dilution of samples mixture with the necessary enrichment media to increase workflow efficiency. This instrument finds its application in mixing in the laboratory or scientific industry.



### KEY FEATURES

- ❖ Speed is adjustable.
- ❖ LCD display indicates the speed and time for BS-RD-Pro.
- ❖ The plate's angle is adjustable from 0° to 90° for BS-RD-Pro.
- ❖ Material of the plate: Stainless steel with oxidation coating (BS-RD-E), Stainless Steel (BS-RD-Pro).
- ❖ Material of the clip: Plastic (BS-RD-E), Stainless Steel (BS-RD-Pro).
- ❖ Three clip sizes are selectable, suit for the tubes: 1.5ml, 15ml, 50ml.



### ACCESSORIES



Cat. No.: 18900160  
Type 1 plate  
60\*1.5ml centrifuge tube



Cat. No.: 18900161(Optional)  
Type 2 plate  
16\*15ml centrifuge tube



Stainless steel clip  
Cat No.: 18900196, for 1.5ml centrifuge tube  
Cat No.: 18900197, for 15ml centrifuge tube  
Cat No.: 18900198, for 50ml centrifuge tube




Plastic clip

Cat No.: 18900327, for 1.5ml centrifuge tube

Cat No.: 18900328, for 15ml centrifuge tube

Cat No.: 18900329, for 50ml centrifuge tube



Cat. No.: 18900162

Type 3 plate

8\*50ml centrifuge tube



Cat. No.: 18900140

Plate support rods

4pcs in a group, use with plates

## Technical Specifications:

Model	BS-RD-E	BS-RD-Pro
Speed Range	0~80rpm	10~70rpm
Inclining Angle Range	/	0°~90°
Timing Range	/	1~1199min
Operation Mode	Continuous	Continuous or Timing
Display	/	LCD
Motor Type	DC motor	
Power Consumption	20W	40W
Power Supply	AC110/220V±10%, 50/60Hz	
Standard Accessories	1 pc Plate, 3 bags Plastic clip (12 pcs/ bag)	1 pc Plate, 1 bag SS clip
External Size (W*D*H)mm	220*300*310	210*280*300
Package Size (W*D*H)mm	350*450*420	350*400*400
Gross Weight(kg)	5	7

### BIOGENIX SYSTEMS PRIVATE LIMITED

S-10, Pankaj Central Plaza, Pocket-5, Plot no-5,

Sector-12, Dwarka, New Delhi-110075, INDIA

Email: [sales@biogenixsystems.com](mailto:sales@biogenixsystems.com) |

Phone: +91-11-43738899 | 9599199767-69 | 9911717057

Website: [www.biogenixsystems.com](http://www.biogenixsystems.com)

Dealer Details

[info@biogenixsystems.com](mailto:info@biogenixsystems.com)

[www.biogenixsystems.com](http://www.biogenixsystems.com)