

Vertical Rotating Mixture Brand - BioGenix®

PRODUCT DESCRIPTION

BIOGENIX® Rotating Mixture is a General purpose instrument which is designed for the dilution of samples mixture with the necessary enrichment media to increase workflow efficiency. This instrument finds its application in mixing in the laboratory or scientific industry.



KEY FEATURES

- ❖ Speed is adjustable
- ❖ Speed range from 0-80 RPM.
- ❖ LCD Display indicates the speed and time for BS-MX-PRO
- ❖ Material of the tube holder: Stainless steel with oxidation coating (BS-MX), Stainless Steel (BS-MX-Pro).
- ❖ Material of the clip: Plastic (BS-MX), stainless Steel (BS-MX-Pro).
- ❖ Three clip sizes are selectable, suit for the tubes: 1.5 ml, 15 ml, 50 ml.
- ❖ A wide range of accessories available.



ACCESSORIES



Cat. No.: 18900142
For 48*1.5ml centrifuge
tube
Held horizontally



Cat. No.: 18900143
(Optional)
For 24*15ml centrifuge
tube
Held horizontally



Cat. No.: 18900144
For 24*50ml centrifuge
tube
Held horizontally



Cat. No.: 18900145
For 32*1.5ml centrifuge tube
Held vertically



Cat. No.: 18900146
(Optional)
For 16*15ml centrifuge tube
Held vertically



Cat. No.: 18900147
For 16*50ml centrifuge tube
Held vertically

Technical Specifications:

Model	BS-MX	BS-MX-Pro (Digital)
Speed Range	0~80rpm	10~70rpm
Timing Range	/	0~1199min
Operation Mode	Continuous	Continuous or Timing
Display	/	LCD
Motor Type	DC motor	
Power Consumption	20W	40W
Power Supply	AC110/220V±10%, 50/60Hz	
Standard Accessories	1 pc Tube holder, 3 bags Plastic clip (12 pcs/ bag)	1pc Tube holder, 1 bag SS clip
External Size (W*D*H)	530*150*190mm	510*220*260mm
Package Size (W*D*H)	800*350*420mm	800*350*420mm
Gross Weight	10kg	11kg

Dealer Details

BIOGENIX SYSTEMS PRIVATE LIMITED

S-10, Pankaj Central Plaza, Pocket-5, Plot no-5,

Sector-12, Dwarka, New Delhi-110075, INDIA

Email: sales@biogenixsystems.com |

Phone: +91-11-43738899 | 9599199767-69 | 9911717057

Website: www.biogenixsystems.com

info@biogenixsystems.com

www.biogenixsystems.com