

Moisture Analyzer

Brand - **BIOGENIX®**

Model – BSMB60. Automation Grade: Automatic

Product Description

Keeping in the mind variegated demands of the clients, we are engaged in offering **BioGenix Moisture Balance**. Our offered balance is especially designed and developed by our vendors for its application in moisture determination and analysis for research as well as quality control purpose. It offers time, auto and manual detection system.



Key Features

- ✓ Halogen light heating
- ✓ Stainless steel chamber
- ✓ Temperature and time can be set
- ✓ The percentage of moisture content
- ✓ The percentage of dry residual
- ✓ LCD with backlight
- ✓ Store historical setting



Soft Touch Keypad

Highlights

- Precision design
- Outstanding performance
- Smooth finish



Moisture analyzer

Application

- Widely used industrial, agriculture, laboratory applications.
- Some typical applications are: tobacco, paper, food, tea, corn, rice, chemical, raw materials and textile materials etc

SPECIFICATION

Model Name/Number	BSMB60 and BSMB110
Automation Grade	Automatic
Usage/Application	Food, Pharma, Industrial Etc.
Brand	BioGenix®
Display Type	Digital
Readability	0.001g or 0.01 percent
Maximum Capacity	60g and 110g
Pan Size	90mm
Size/Dimension	240 x365 x180
Temperature Range	40 deg to 199 deg.
Resolution	0.01 percent
Data Memory Function Type	15
Heating Technology	Halogen Lamp
Weight	4.82 kg
Communication	RS232
Minimum Order Quantity	1

Advantages

- ☐ Quick measurement
- ☐ Portable
- ☐ The halogen moisture analyzer does not require any calculation, and
- ☐ its sampling is simple; therefore, the possibility of errors decreases.
- ☐ The heating rate is fast; therefore, the halogen moisture analyzer can
- ☐ quickly reach its maximum power
- ☐ The equipment is fully automated
- ☐ Needs a very small space

Dealer Details



BIOGENIX SYSTEMS PRIVATE LIMITED

S-10, Pankaj Central Plaza, Pocket-5, Plot no-5,
Sector-12, Dwarka, New Delhi-110075, INDIA

Email: sales@biogenixsystems.com |

Phone: +91-11-43738899 | 9911717057

Website: www.biogenixsystems.com

info@biogenixsystems.com
www.biogenixsystems.com